



## Roll-On Power to AMP Up Your Life!

**AMP's** unique, patented design helps take you where you want to be - even when you don't have the energy or endurance. Keep the lightweight benefits of your manual wheelchair and only add on power when and where you need it. When you are fatigued or find the terrain challenging – add that extra boost of power!

This add-on device is powered by a lightweight, powerful Lithium Ion Phosphate battery (LiFeOP4). Batteries with twice the power, more cycles and faster charging than typical lead acid batteries. It will keep you going farther for extended time – up to 11 miles on one charge!

**AMP** is easy – easy to dock on and release off. Best of all, its intuitive to use!

While seated in your wheelchair, just roll on or off **AMP** – it docks or undocks with a click. Control your **AMP** with an easy-to-use, instinctive joystick. The joystick can also be set up for attendant driving too.

Let **AMP** give you that extra help when you need it - let **AMP** take you where you want to be!

AMP by Method Mobility...for where you want to be.

# AMP Specifications

- **Available for:** Most rigid, folding or tilt wheelchairs with 24" or 25" wheels.\*
- **Wheelchair Width Compatibility:** 14"-20"
- **Wheelchair Rear Lower Frame Height Compatibility:** Frame must have at least 8.5" ground clearance.
- **Power Module Weight:** Less than 42 lbs. with battery. Battery Weight: 8 lbs.
- **Weight Limit:** 265 lbs.
- **Range:** Up to 11 miles on single charge.
- **Speed:** Up to 4 mph
- **Suitable for:** Use indoors or outdoors on ADA specified surfaces.



1. Roll on, lock wheels, you are docked!
2. Plug in and turn on!
3. Go! Go! Go!



Attendant Mount or  
Swing-Away Joystick



**Method Mobility**  
2562 N. Fordham Avenue  
Fresno, CA 93727 USA  
Telephone: (888) 88AMP88 (882-6788)  
Email: [CustomerCare@methodmobility.com](mailto:CustomerCare@methodmobility.com)