

Mark For: _____
 Dealer/VA Acct. #: _____
 Dealer/VA Name: _____
 Dealer/VA Contact Phone: _____
 Dealer/VA City: _____
 Order Confirmation Email: _____

Quote Order

P.O. # _____

Date: _____

Ship to: _____

Address: _____

City, State, Zip: _____






ORDER FORM

WHEELCHAIR SPECIFICATIONS

End User Wheelchair Type: Folding Rigid Tilt Brand _____ Model _____

End User Wheelchair Specs^{1, 2}:



¹Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". ²To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. ³At rear of chair inside distance between tires at hub height.

AMP OPTIONS

- \$7,995.00
- Joystick Input (Armrest Mount) Included
- Must choose one:
 - Armrest Swing-Away Mount \$350.00
 - Attendant Control Joystick Mount (with AMP Joystick Input) \$440.00
 - Side: Right Left (must choose one)
- Joystick Handles
 - Standard Handle Included
 - Bodypoint® 3" U-Shaped Handle \$135.00
 - Bodypoint® Grooved Mushroom Handle \$135.00
- Lithium Iron Phosphate LiFeP04 Battery and Charger Included

AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

Check One

AMP SIZE MATRIX		
FOLDING/RIGID 0/3 CAMBER		
CHAIR SEAT WIDTH		AMP SIZE (growth)
14"		SMALL (14"-16")
15"		SMALL (14"-16")
16"		SMALL (14"-16")
16"		MEDIUM (16"-18")
17"		MEDIUM (16"-18")
18"		MEDIUM (16"-18")
18"		LARGE (18"-20")
19"		LARGE (18"-20")
20"		LARGE (no growth)

BATTERY OPTIONS

- Back Up Lithium Iron Phosphate LiFeP04 Battery \$950.00
- Spare Battery Charger \$225.00

MILITARY BRANCH LOGO (VA ONLY)

- Army Navy Air Force Coast Guard Marines No Charge

Check One

AMP SIZE MATRIX		
FOLDING/RIGID >5 CAMBER TILT 0/3 CAMBER		
CHAIR SEAT WIDTH		AMP SIZE (growth)
14"		SMALL (14"-15")
15"		MEDIUM (15"-17")
16"		MEDIUM (15"-17")
17"		LARGE (17"-19")
18"		LARGE (17"-19")
19"		LARGE (no growth)
20"	N/A	

NOTE: All Specifications are subject to change.
Always refer to an Authorized Method Mobility supplier.



888-88AMP88 (882-6788)
 www.methodmobility.com

For quotes and orders:
 CustomerCare@methodmobility.com

MM-AMPOF 6.23.6