



Taking You Where You Want to Be

AMP Power Add-On



Easy Roll-On



Docked

Easy Roll-Off



Undocked

Retractable Joystick



AMP Specifications

- **Available for:** Most rigid, folding or tilt wheelchairs with 22", 24", 25" or 26" wheels. Some restrictions apply.
- **Wheelchair Width Compatibility:** 14"-20"
- **Power Module Weight:** Less than 42 lbs. with battery. **Battery Weight:** 8 lbs.
- **PG Controller with Joystick Mounts.**
- **Weight Limit:** 300 lbs.
- **Unique LifePO4 Lithium Iron Phosphate Battery.**
- **Range:** Up to 11 miles on single charge.
- **Speed:** Up to 4 mph

Your Method Mobility Sales Representative



2562 N. Fordham Avenue
Fresno, CA 93727 USA
Telephone: (888) 88AMP88 (882-6788)
Email: CustomerCare@methodmobility.com