

# AMP JOYSTICK

Mark For: \_\_\_\_\_  
 Acct #: \_\_\_\_\_  
 Bill To Name: \_\_\_\_\_  
 Billing Address: \_\_\_\_\_  
 City, State: \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Contact : \_\_\_\_\_

QUOTE ☐ ORDER ☐ #: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 PO#: \_\_\_\_\_  
 Ship to: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, State: \_\_\_\_\_  
 Zip: \_\_\_\_\_



## END USER WHEELCHAIR SPECIFICATIONS

Wheelchair Type: Folding ☐ Rigid ☐ Tilt ☐ Brand: \_\_\_\_\_ Model: \_\_\_\_\_  
 1. Front Seat Height in inches \_\_\_\_\_  
 2. Rear Seat Height in inches \_\_\_\_\_  
 3. Inside Tire Distance in inches \_\_\_\_\_

<sup>1</sup> Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". <sup>2</sup> To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. <sup>3</sup> At bottom of chair inside to inside distance between tires.

## JOYSTICK AMP & OPTIONS

| SKU                             | Description   | Retail  |
|---------------------------------|---|---------|
| <input type="checkbox"/> 180AJS | AMPjs Power Module w/ Battery, Charger and Joystick | \$7,595 |

(Must choose one from the size matrix to the right:)

- ☐ 180AJSS AMP Small Base 14" - 16" Wide  
☐ 180AJSM AMP Medium Base 16" - 18" Wide  
☐ 180AJSL AMP Large Base 18" - 20" Wide

☐ 180WEXT Width Extension Kit (Adds 2" to Wheel Width Fit) \$275

Joystick Mount Side (Must choose one)

- ☐ 180JIR Joystick Input Right  
☐ 180JIL Joystick Input Left

Swing Away/Attendant Control (Must choose at least one)

- ☐ 180SAD Swing Away Mount Dual Post \$350  
☐ 180SAT Swing Away Mount Tubular \$350  
☐ 180ACP Attendant Control Padded \$440  
☐ 180ACF Attendant Control Folding \$440

### Joystick Handles

- ☒ 180SH Standard Handle Included  
☐ 180BPU Bodypoint® 3" U-Shaped Handle \$135  
☐ 180BPGM Bodypoint® Grooved Mushroom Handle \$135

### Battery/Charger

- ☒ 180LPSTC LiFePO4 Battery /Charger (6.5 Hours) Included  
☐ 180USC Upgrade Super Charger (2.5 hours) \$125

### Battery Options

- ☐ 180BULP Back Up LiFePO4 Battery \$995  
☐ 180BUSEC Spare Battery Charger (6.5 Hours) \$225  
☐ 180SSC Spare Super Charger (2.5 Hours) \$335

### Military Branch Logo (VA ONLY)

- ☐ Army NC  
☐ Navy NC  
☐ Air Force NC  
☐ Marines NC  
☐ Coast Guard NC

Sub Totals

## AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

### CHECK ONE BELOW

| FOLDING/RIGID 0/3 CAMBER |  | WHEEL WIDTH RANGE<br> |
|--------------------------|--|-----------------------|
| CHAIR SEAT WIDTH         | AMP SIZE (growth)                          |                       |
| 14"                      | <input type="checkbox"/> SMALL (14"-16")   | 17 1/4"-19 1/2"       |
| 15"                      | <input type="checkbox"/> SMALL (14"-16")   |                       |
| 16"                      | <input type="checkbox"/> SMALL (14"-16")   |                       |
| 16"                      | <input type="checkbox"/> MEDIUM (16"-18")  | 19" - 21 1/4"         |
| 17"                      | <input type="checkbox"/> MEDIUM (16"-18")  |                       |
| 18"                      | <input type="checkbox"/> MEDIUM (16"-18")  |                       |
| 18"                      | <input type="checkbox"/> LARGE (18"-20")   | 20 1/2" - 23"         |
| 19"                      | <input type="checkbox"/> LARGE (18"-20")   |                       |
| 20"                      | <input type="checkbox"/> LARGE (No Growth) |                       |

| FOLDING/RIGID UP TO 5 CAMBER<br>TILT 0 CAMBER |  | WHEEL WIDTH RANGE<br> |
|---|--|-----------------------|
| CHAIR SEAT WIDTH                              | AMP SIZE (growth)                          |                       |
| 14"   | <input type="checkbox"/> SMALL (14"-15")   | 17 1/4"-19 1/2"       |
| 15"   | <input type="checkbox"/> MEDIUM (15"-17")  |                       |
| 16"   | <input type="checkbox"/> MEDIUM (15"-17")  |                       |
| 17"   | <input type="checkbox"/> LARGE (17"-19")   | 20 1/2" - 23"         |
| 18"   | <input type="checkbox"/> LARGE (17"-19")   |                       |
| 19"   | <input type="checkbox"/> LARGE (No Growth) |                       |
| 20"   | N/A  |                       |

For quotes and orders EMAIL to:  
 CustomerCare@methodmobility.com  
 or FAX to: 877-398-4545

Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9

www.methodmobility.com  
 888-88AMP88 (882-6788)

|                                       |        |
|---------------------------------------|--------|
| Shipping                              |        |
| Total w/ Assess., Delivery & Shipping | \$0.00 |
|                                       | \$0.00 |
| Total with Fees(Pay this amount)      | \$0.00 |



MM-AMPOF JS 8.25.1