

# AMP HANDRIM

Mark For: \_\_\_\_\_  
 Acct #: \_\_\_\_\_  
 Bill To Name: \_\_\_\_\_  
 Billing Address: \_\_\_\_\_  
 City, State: \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Contact : \_\_\_\_\_

QUOTE ☐ ORDER ☐ #: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 PO#: \_\_\_\_\_  
 Ship to: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, State: \_\_\_\_\_  
 Zip: \_\_\_\_\_



## END USER WHEELCHAIR SPECIFICATIONS

Wheelchair Type: Folding ☐ Rigid ☐

Brand: \_\_\_\_\_

Model: \_\_\_\_\_

1. Front Seat Height in inches




2. Rear Seat Height in inches




3. Inside Tire Distance in inches




1 Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". 2 To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. 3 At bottom of chair inside to inside distance between tires.

## AMP HANDRIM & OPTIONS

☐ AMP E0986 AMP Handrim Power Module w/ Battery & Charger \$7,995

(Must choose one from the size matrix to the right:)

- ☐ 380AHR5 AMP Small Base 14" - 16" Wide  
☐ 380AHRM AMP Medium Base 16" - 18" Wide  
☐ 380AHLR AMP Large Base 18" - 20" Wide  
☐ 180WEXT Width Extension Kit (Adds 2" to Wheel Width Fit) \$275

### Joystick Input

☐ 380JIO Joystick Input \$595

Joystick Mount Side (Must choose one if ordering Joystick option)

- ☐ 380JIR Joystick Input Right  
☐ 380JIL Joystick Input Left

Swing Away/Attendant Control (Must choose at least one if ordering Joystick option)

- ☐ 380SAD Swing Away Mount Dual Post \$350  
☐ 380SAT Swing Away Mount Tubular \$350

### Joystick Handles

- ☒ 380SH Standard Handle (When ordering Joystick input) Included  
☐ 380BPU Bodypoint® 3" U-Shaped Handle \$135  
☐ 380BPGM Bodypoint® Grooved Mushroom Handle \$135

### Battery/ Charger

- ☒ 380LPSTC LiFePO4 Battery/ Charger (6.5 Hours) Included  
☐ 380USC Upgrade Super Charger (2.5 hours) \$125

### Battery Options

- ☐ 380BULP Back Up LiFePO4 Battery \$995  
☐ 380BUSE Spare Battery Charger (6.5 Hours) \$225  
☐ 380SSC Spare Super Charger (2.5 Hours) \$335

### Military Branch Logo (VA ONLY)

- ☐ Army NC  
☐ Navy NC  
☐ Air Force NC  
☐ Marines NC  
☐ Coast Guard NC

Sub Totals

## AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

### CHECK ONE BELOW

FOLDING/RIGID 0/3 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	<input type="checkbox"/> SMALL (14"-16")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> MEDIUM (16"-18")	19" - 21 1/4"
17"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> LARGE (18"-20")	20 1/2" - 23"
19"	<input type="checkbox"/> LARGE (18"-20")	
20"	<input type="checkbox"/> LARGE (No Growth)	

FOLDING/RIGID UP TO 5 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	<input type="checkbox"/> SMALL (14"-15")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> MEDIUM (15"-17")	
16"	<input type="checkbox"/> MEDIUM (15"-17")	
17"	<input type="checkbox"/> LARGE (17"-19")	20 1/2" - 23"
18"	<input type="checkbox"/> LARGE (17"-19")	
19"	<input type="checkbox"/> LARGE (No Growth)	
20"	N/A	

For quotes and orders EMAIL to:  
 CustomerCare@methodmobility.com  
 or FAX to: 877-398-4545

Shipping  
 Total w/Assess.,Delivery & Shipping  
 Credit Card Fee  
 Total with Fees (Pay this amount)



Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9

www.methodmobility.com  
 888-88AMP88 (882-6788)

MM-AMPOF 8.25.1