

AMP JOYSTICK

Mark For: _____
 Acct #: _____
 Bill To Name: _____
 Billing Address: _____
 City, State: _____
 Confirmation Email: _____
 Contact : _____

QUOTE ☐ ORDER ☐ #: _____
 Date: _____
 PO#: _____
 Ship to: _____
 Address: _____
 City, State: _____
 Zip: _____



END USER WHEELCHAIR SPECIFICATIONS

Wheelchair Type: Folding ☐ Rigid ☐ Tilt ☐ Brand: _____ Model: _____
 1. Front Seat Height in inches _____
 2. Rear Seat Height in inches _____
 3. Inside Tire Distance in inches _____

¹ Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". ² To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. ³ At bottom of chair inside to inside distance between tires.

JOYSTICK AMP & OPTIONS

SKU	HPCS	Description	Retail
<input type="checkbox"/> 180AJS	E0986	AMPjs Power Module w/ Battery, Charger and Joystick	\$7,595

(Must choose one from the size matrix to the right:)

- ☐ 180AJSS AMP Small Base 14" - 16" Wide
☐ 180AJSM AMP Medium Base 16" - 18" Wide
☐ 180AJSL AMP Large Base 18" - 20" Wide

☐ 180WEXT Width Extension Kit (Adds 2" to Wheel Width Fit) \$275

Joystick Mount Side (Must choose one)

- ☐ 180JIR Joystick Input Right
☐ 180JIL Joystick Input Left

Swing Away/Attendant Control (Must choose at least one)

- ☐ 180SAD Swing Away Mount Dual Post \$350
☐ 180SAT Swing Away Mount Tubular \$350
☐ 180ACP Attendant Control Padded \$440
☐ 180ACF Attendant Control Folding \$440

Joystick Handles

- ☒ 180SH Standard Handle Included
☐ 180BPU Bodypoint® 3" U-Shaped Handle \$135
☐ 180BPGM Bodypoint® Grooved Mushroom Handle \$135

Battery/Charger

- ☒ 180LPSTC LiFePO4 Battery /Charger (6.5 Hours) Included
☐ 180USC Upgrade Super Charger (2.5 hours) \$125

Battery Options

- ☐ 180BULP Back Up LiFePO4 Battery \$995
☐ 180BUSEC Spare Battery Charger (6.5 Hours) \$225
☐ 180SSC Spare Super Charger (2.5 Hours) \$335

Military Branch Logo (VA ONLY)

- ☐ 180ABL Army NC
☐ 180NBL Navy NC
☐ 180AFBL Air Force NC
☐ 180MBL Marines NC
☐ 180CGBL Coast Guard NC

Sub Totals

For quotes and orders EMAIL to:
 CustomerCare@methodmobility.com
 or FAX to: 877-398-4545

Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9

www.methodmobility.com
 888-88AMP88 (882-6788)

Shipping	
Total w/ Assess.,Delivery & Shipping	\$0.00
	\$0.00
Total with Fees(Pay this amount)	\$0.00

AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

CHECK ONE BELOW

FOLDING/RIGID 0/3 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	
14"	<input type="checkbox"/> SMALL (14"-16")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> MEDIUM (16"-18")	19" - 21 1/4"
17"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> LARGE (18"-20")	20 1/2" - 23"
19"	<input type="checkbox"/> LARGE (18"-20")	
20"	<input type="checkbox"/> LARGE (No Growth)	

FOLDING/RIGID UP TO 5 CAMBER TILT 0 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	
14"	<input type="checkbox"/> SMALL (14"-15")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> MEDIUM (15"-17")	
16"	<input type="checkbox"/> MEDIUM (15"-17")	
17"	<input type="checkbox"/> LARGE (17"-19")	20 1/2" - 23"
18"	<input type="checkbox"/> LARGE (17"-19")	
19"	<input type="checkbox"/> LARGE (No Growth)	
20"	N/A	



MM-AMPOF JS 11.5.1