

# AMP HANDRIM

Mark For: _____ Acct #: _____ Bill To Name: _____ Billing Address: _____ City, State: _____ Confirmation Email: _____ Contact : _____	QUOTE <input type="checkbox"/> ORDER <input type="checkbox"/> #: _____ Date: _____ PO#: _____ Ship to: _____ Address: _____ City, State: _____ Zip: _____	
---	---	--

**END USER WHEELCHAIR SPECIFICATIONS**

Wheelchair Type: Folding  Rigid  Brand: \_\_\_\_\_ Model: \_\_\_\_\_

1. Front Seat Height in inches \_\_\_\_\_

2. Rear Seat Height in inches \_\_\_\_\_

3. Inside Tire Distance in inches \_\_\_\_\_

1 Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". 2 To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. 3 At bottom of chair inside to inside distance between tires.

AMP HANDRIM & OPTIONS			
SKU	HCPC	Description	Retail
<input type="checkbox"/>	AMP E0986	AMP Handrim Power Module w/ Battery & Charger	\$7,995
<b>(Must choose one from the size matrix to the right:)</b>			
<input type="checkbox"/>	380AHRS	AMP Small Base 14" - 16" Wide	
<input type="checkbox"/>	380AHRM	AMP Medium Base 16" - 18" Wide	
<input type="checkbox"/>	380AHL	AMP Large Base 18" - 20" Wide	
<input type="checkbox"/>	180WEXT	Width Extension Kit (Adds 2" to Wheel Width Fit)	\$275
<b>Joystick Input</b>			
<input type="checkbox"/>	380JIO	Joystick Input	\$595
Joystick Mount Side <b>(Must choose one if ordering Joystick option)</b>			
<input type="checkbox"/>	380JIR	Joystick Input Right	
<input type="checkbox"/>	380JIL	Joystick Input Left	
Swing Away/Attendant Control <b>(Must choose at least one if ordering Joystick option)</b>			
<input type="checkbox"/>	380SAD	Swing Away Mount Dual Post	\$350
<input type="checkbox"/>	380SAT	Swing Away Mount Tubular	\$350
<b>Joystick Handles</b>			
<input checked="" type="checkbox"/>	380SH	Standard Handle (When ordering Joystick input)	Included
<input type="checkbox"/>	380BPU	Bodypoint® 3" U-Shaped Handle	\$135
<input type="checkbox"/>	380BPGM	Bodypoint® Grooved Mushroom Handle	\$135
<b>Battery/ Charger</b>			
<input checked="" type="checkbox"/>	380LPSTC	LiFePO4 Battery/ Charger (6.5 Hours)	Included
<input type="checkbox"/>	380USC	Upgrade Super Charger (2.5 hours)	\$125
<b>Battery Options</b>			
<input type="checkbox"/>	380BULP	Back Up LiFePO4 Battery	\$995
<input type="checkbox"/>	380SBC	Spare Battery Charger (6.5 Hours)	\$225
<input type="checkbox"/>	380SSC	Spare Super Charger (2.5 Hours)	\$335
<b>Military Branch Logo (VA ONLY)</b>			
<input type="checkbox"/>	380ABL	Army	NC
<input type="checkbox"/>	380NBL	Navy	NC
<input type="checkbox"/>	380AFBL	Air Force	NC
<input type="checkbox"/>	380MBL	Marines	NC
<input type="checkbox"/>	380CGBL	Coast Guard	NC
<b>Sub Totals</b>			

AMP SIZE MATRIX
NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

CHECK ONE BELOW		
FOLDING/RIGID 0/3 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	<input type="checkbox"/> SMALL (14"-16")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> MEDIUM (16"-18")	19" - 21 1/4"
17"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> LARGE (18"-20")	20 1/2" - 23"
19"	<input type="checkbox"/> LARGE (18"-20")	
20"	<input type="checkbox"/> LARGE (No Growth)	

FOLDING/RIGID UP TO 5 CAMBER		
CHAIR SEAT WIDTH	AMP SIZE (growth)	WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	<input type="checkbox"/> SMALL (14"-15")	17 1/4"-19 1/2"
15"	<input type="checkbox"/> MEDIUM (15"-17")	
16"	<input type="checkbox"/> MEDIUM (15"-17")	
17"	<input type="checkbox"/> LARGE (17"-19")	19" - 21 1/4"
18"	<input type="checkbox"/> LARGE (17"-19")	
19"	<input type="checkbox"/> LARGE (No Growth*)	
20"	N/A	20 1/2" - 23"

For quotes and orders EMAIL to:  
**CustomerCare@methodmobility.com**  
 or FAX to: 877-398-4545

Shipping	
Total w/Assess., Delivery & Shipping	
Credit Card Fee	
<b>Total with Fees (Pay this amount)</b>	



**Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9**

www.methodmobility.com  
 888-88AMP88 (882-6788)