

AMP JOYSTICK



Mark For: _____
 Acct #: _____
 Bill To Name: _____
 Billing Address: _____
 City, State: _____
 Confirmation Email: _____
 Contact : _____

QUOTE ORDER #: _____
 Date: _____
 PO#: _____
 Ship to: _____
 Address: _____
 City, State: _____
 Zip: _____

END USER WHEELCHAIR SPECIFICATIONS

Wheelchair Type: Folding Rigid Tilt Brand: _____ Model: _____

1. Front Seat Height in inches _____ 2. Rear Seat Height in inches _____ 3. Inside Tire Distance in inches _____

¹ Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". ² To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. ³ At bottom of chair inside to inside distance between tires.

JOYSTICK AMP & OPTIONS

SKU HCPC	Description	Retail
<input type="checkbox"/> 180AJS E0986	AMPjs Power Module w/ Battery, Charger and Joystick	\$7,595
(Must choose one from the size matrix to the right:)		
<input type="checkbox"/> 180AJSS	AMP Small Base 14" - 16" Wide	
<input type="checkbox"/> 180AJSM	AMP Medium Base 16" - 18" Wide	
<input type="checkbox"/> 180AJSL	AMP Large Base 18" - 20" Wide	
<input type="checkbox"/> 180WEXT	Width Extension Kit (Adds 2" to Wheel Width Fit)	\$275
Joystick Mount Side (Must choose one)		
<input type="checkbox"/> 180JIR	Joystick Input Right	
<input type="checkbox"/> 180JIL	Joystick Input Left	
Swing Away/Attendant Control Mount (Must choose at least one)		
<input type="checkbox"/> 180SAD	Swing Away Mount Dual Post	\$350
<input type="checkbox"/> 180SAT	Swing Away Mount Tubular	\$350
<input type="checkbox"/> 180ACP	Attendant Control Mount/Padded Push Handle	\$440
<input type="checkbox"/> 180ACF	Attendant Control Mount /Folding Push Handle	\$440
Additional Joystick when utilizing two inputs above		
<input type="checkbox"/> 180ADJS	Additional Joystick control	\$650
Joystick Handles		
<input checked="" type="checkbox"/> 180SH	Standard Handle	Included
<input type="checkbox"/> 180BPU	Bodypoint® 3" U-Shaped Handle	\$135
<input type="checkbox"/> 180BPGM	Bodypoint® Grooved Mushroom Handle	\$135
Battery/Charger		
<input checked="" type="checkbox"/> 180LPSTC	LiFePO4 Battery /Charger (6.5 Hours)	Included
<input type="checkbox"/> 180USC	Upgrade Super Charger (2.5 hours)	\$125
Battery Options		
<input type="checkbox"/> 180BULP	Back Up LiFePO4 Battery	\$995
<input type="checkbox"/> 180SBC	Spare Battery Charger (6.5 Hours)	\$225
<input type="checkbox"/> 180SSC	Spare Super Charger (2.5 Hours)	\$335
Military Branch Logo (VA ONLY)		
<input type="checkbox"/> 180ABL	Army	NC
<input type="checkbox"/> 180NBL	Navy	NC
<input type="checkbox"/> 180AFBL	Air Force	NC
<input type="checkbox"/> 180MBL	Marines	NC
<input type="checkbox"/> 180CGBL	Coast Guard	NC

Sub Totals

For quotes and orders EMAIL to:
CustomerCare@methodmobility.com
 or FAX to: 877-398-4545

Shipping	
Total w/ Assess., Delivery & Shipping	\$0.00
Total with Fees (Pay this amount)	\$0.00

Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9

www.methodmobility.com

AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

CHECK ONE BELOW

FOLDING/RIGID 0/3 CAMBER		WHEEL WIDTH RANGE INSIDE TO INSIDE OF TIRES
CHAIR SEAT WIDTH	AMP SIZE (growth)	
14"	<input type="checkbox"/> SMALL (14"-16")	17 1/4" - 19 1/2"
15"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> SMALL (14"-16")	
16"	<input type="checkbox"/> MEDIUM (16"-18")	19" - 21 1/4"
17"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> MEDIUM (16"-18")	
18"	<input type="checkbox"/> LARGE (18"-20")	20 1/2" - 23"
19"	<input type="checkbox"/> LARGE (18"-20")	
20"	<input type="checkbox"/> LARGE (No Growth)	

FOLDING/RIGID UP TO 5 CAMBER TILT 0 CAMBER		WHEEL WIDTH RANGE INSIDE TO INSIDE OF TIRES
CHAIR SEAT WIDTH	AMP SIZE (growth)	
14"	<input type="checkbox"/> SMALL (14"-15")	17 1/4" - 19 1/2"
15"	<input type="checkbox"/> MEDIUM (15"-17")	
16"	<input type="checkbox"/> MEDIUM (15"-17")	
17"	<input type="checkbox"/> LARGE (17"-19")	20 1/2" - 23"
18"	<input type="checkbox"/> LARGE (17"-19")	
19"	<input type="checkbox"/> LARGE (No Growth)	
20"	N/A	



MM-AMPOF JS 2.26.1